




JOFEL Antiseptic Cleansing Lotion

✓ **Targets markets:**
food processing, kitchens and institutional catering

- Broad spectrum of action
- Dermatologically tested
- Non-coloured, unscented
- Enriched with glycerine
- Complies with regulation 528/2012



TECHNICAL INFORMATION

Substances actives



- Chlorhexidine digluconate CAS n° 18472-51-0 (0.1% m/m)
- Didecylidimethylammonium chloride CAS 7173-51-5 (0.85% m/m)

Caractéristiques



Appearance : translucid creamy liquid
Fragrance : none
Colour : colourness
pH : 6,0 – 8,0

Innocuousness



- Dermatologically tested (non-irritating)
- Designed for frequent washing
- Clinically proven to help keep skin healthy

Conservation



- The product will keep for 24 months in its original packaging under normal storage conditions. CRAINT LE GEL, STOCKER ENTRE 5 ET 35°C
- Protect from frost, store between 5 and 35°C.

Instructions for use



- Run your hands under water.
- Place a dose of product on the palm of your hand.
- Soap your hands, interdigital spaces and around your nails for at least 30 seconds.
- Rinse hands thoroughly, then dry with a single-use towel.

Réglementation



- TP1 biocidal product
- For professional use only

Environment



- Strict requirements met with regard to harmful substances
- Reduction of packaging waste.

Precautions for use



- Do not swallow. Do not mix with other products.
- Use biocides with care.
- Read the label and product information before use.
- The packaging must be disposed of as hazardous waste, or non-hazardous if rinsed, under the sole responsibility of the waste holder.
- Do not discharge the pure product into the environment.
- Do not reuse empty packaging.
- For more information, refer to the product's SDS
- In an emergency, contact the nearest poison control centre.

MICROBIOLOGICAL PROPERTIES

Active on	Standards	Contact time
Bacteria	EN 1276* (v. August 2019) EN 1499 (v. 2013), + active against Escherichia coli K12	60s 30s 30s
Yeast	EN 1650* (v. August 2019), active against Candida albicans	60s
Virus	EN 14476* (v. July 2019), active against vaccinia virus Spectrum of activity covering enveloped viruses including: Orthopoxvirus (including the virus responsible for mmonkeypox), Herpes, Influenza A (H1N1), Hepatitis B, Hepatitis C, HIV...	60s

*Normes réalisées en conditions de saleté / Standards achieved in dirty conditions

CONDITIONING	DISTRIBUTION SYSTEM	PRODUCT CODE
Recharge AITANA 1000mL x 6	Distributeur AITANA MIX (manuel)	AT54906F
Recharge TIMELESS 1000mL x 6	Distributeur TIMELESS (manuel)	AT20906F
Recharge OPTIQUE 1000mL x 6	Distributeur AITANA OPTIQUE (automatique)	AT18906F
Recharge OPTIQUE 1000mL x 6	Distributeur TIMELESS OPTIQUE (automatique)	AT18906F

PALLET MEASUREMENTS



Palette

80cm x 120 cm x 170cm



Total

78 boxes