




JO FEL Hydroalcoholic Gel

✓ **Target Markets :**
For use in medical, paramedical, veterinary, food processing, food service, childcare, and other similar environments.

- Biocidal product intended for human hygiene (TP1).
- For rapid antisepsis of the skin and hands.



TECHNICAL INFORMATION

Characteristics

Appearance: Gel
Smell: Alcohol
Color: Colorless
Density: 0,891

Instructions for use

- The HYDROALCOHOLIC GEL disinfects in one step: hand hygiene is ensured
- Apply the gel to dry, clean hands
- Do not rinse, do not use a towel
- Use regularly, as hands are a recognized vehicle for microbial contamination
- After use, the HYDROALCOHOLIC GEL leaves a pleasant feeling of freshness with no smell

Safety Instructions

- HYDROALCOHOLIC GEL is intended for professional use only.
- Do not ingest.
- Avoid contact with eyes. In case of accidental splashes, rinse thoroughly with water.

MICROBIOLOGICAL PROPERTIES

Bactericidal activity	- UNE EN 1276 (60 s) - UNE EN 1500 (60 s) - UNE EN 1040, bactericidal against <i>Pseudomonas aeruginosa</i> and <i>Staphylococcus aureus</i> strains.
Yeasticidal activity	- UNE EN 1650 (60 s)
Virucidal activity	- NF EN 14476 (120 s), virucidal activity against Poliovirus Type 1 and Adenovirus Type 5.

INGREDIENTS

EU INCI name of the ingredient*

AQUA
ALCOHOL
ISOPROPYL ALCOHOL
BENZALKONIUM CHLORIDE
PROPYLENE GLYCOL,
TRIETHANOLAMINE
CYCLOHEXANE
ALCOHOL
C11-14
ETOXILADO

*Nomenclature internationale d'ingrédients cosmétiques

CONDITIONING	DISTRIBUTION SYSTEM	PRODUCT CODE
AITANA 1000mL x 6 Recharge	AITANA MIX Dispenser (manuel)	AT86006F
TIMELESS 1000mL x 6 Recharge	TIMELESS Dispenser (manuel)	AT93006F
OPTICAL 1000mL x 6 Recharge	AITANA OPTICAL Dispenser (automatic)	AT83006F
OPTICAL 1000mL x 6 Recharge	TIMELESS OPTICAL Dispenser (automatic)	AT83006F

PALLET MEASUREMENTS



Palett

80cm x 120 cm x 170cm



Total

78 boxes